

CMTA

Comprehensive Molecular Tumor Analysis

Laboratory medical accreditation in accordance with
DIN EN ISO 15189:2024

The personalized cancer diagnostics service from
Alacris Theranostics

Alacris CMTA® Highlights:

- Concise, interpreted report intended to inform the oncologist for navigating therapy choice
- Integrated deep profiling from whole exome and bulk RNAseq for optimized target identification
- Tailored solution to align with the patient's individual journey

INFORMATION PACKAGE CONTENT

Service Description	03
Genetic Investigation Requisition	04
Declaration of Consent	05
Patient Information Sheet	06
Sample Requirements & Shipping Guidelines	07
Tumor Sample Info Sheet	09
Checklist	10



Best-in-class solution for precision oncology

We harness the full power of whole exome/genome and bulk RNAseq to broaden the range of therapy recommendations.

Every patient is unique, every tumor is different. To capture this individuality, the CMTA technology integrates DNA and RNA information in a comprehensive multi-modal tumor profiling to unleash more actionable insights.

The Alacris CMTA® translates our long-standing expertise in cancer genomics into cancer diagnostic solutions. Our end-to-end CMTA service uses state-of-the-art NGS technologies coupled with proprietary bioinformatics analysis, designed to suit every tumor type.

Our DIN EN ISO 15189:2024 medically accredited laboratory carries out the CMTA test within its accredited scope. Note: the list of genes included in the accreditation scope is available on our website (www.alacris.de).

Comprehensive, unbiased, multi-modal insights

- Mutations, TMB
- Gene fusions
- Expressed variants
- Oncogenic splice variants
- HRD
- Prognostic gene signatures
- Copy number variants (CNVs)
- Immune microenvironment
- Loss of heterozygosity (LOH)
- Deregulated gene expression

Interpreted key features with matched therapeutic options

Why choose the Alacris CMTA®?


Accurate, sensitive, comprehensive test in a single process

Tumor-agnostic test enabling personalized care for all cancers

- Customized to fit the patient's journey
- Flexible analysis for specific needs
- Value-based diagnostics, outstanding data quality
- Data protection in compliance with GDPR regulations
- Competitive pricing

Alacris-CMTA® Workflow


Tumor-agnostic comprehensive profiling




Tissue Sample



NGS Deep Sequencing



Bioinformatics Analysis & Patient Report Generation



From sample to final report

Contact us:

cmtanalysis@alacris.de | www.alacris.de
 Max-Planck-Str.3 | D-12489 Berlin | Germany
 Tel +49 30 8431225 - 10 / - 30 | Fax +49 30 8431225 - 40

Service Description:

Comprehensive Molecular Tumor Analysis

A. Laboratory Analysis Components

1. Isolation of genomic DNA and total RNA from the tumor sample
2. Isolation of genomic DNA from blood/ control tissue
3. Preparation of sequencing libraries:
 - Whole-exome libraries from tumor and blood DNA
 - RNA libraries (ribo-depleted, strand specific) from tumor total RNA
4. Next Generation Sequencing (NGS):
 - Whole exome seq: Paired-end sequencing 100 bp, 18-24 Gb (200-250x mean coverage) of filtered data per library
 - Bulk RNAseq: Paired-end sequencing, 100 bp, 18-20 Gb of filtered data per library

B. Bioinformatics Analysis

1. Mutational burden (TMB) (number of mutations per Megabase), graphical representation of the TMB of the patient in comparison to a cohort from a similar tumor type
2. Somatic variants: mutations (SNVs and Indels) and expressed mutations
3. Mutational signatures
4. Chromosomal copy number variants (CNV: amplifications, deletions, loss of heterozygosity-LOH)
5. Gene fusions
6. Oncogenic splice variants
7. Integrative gene expression profiling: (a) dysregulated genes of clinical interest (b) Expression levels of key biomarkers (c) gene expression in CNV regions, (d) gene signatures analysis on demand (e.g. signaling pathways, drug sensitivity)
8. Estimation of the tumor immune microenvironment

C. Patient Report

Generation of a comprehensive patient report summarizing the key results of the integrative molecular tumor analysis, including interpreted clinically relevant annotations and matched potential treatments.

Note:

If the quantity or quality of the samples (e.g. < 30% tumor content) is insufficient, the analysis cannot be performed appropriately. This only becomes apparent during quality controls within the work process. In this case, only the services provided up to that point will be invoiced.

Terms and conditions

Alacris AGB are valid for this offer and are available on company's webpage:

www.alacris.de/terms-and-conditions

These analyses are intended as information for the treating physician only and - if requested by the patient from their physician - for the patient, or in the case of express and prior written consent, also for certain third parties appointed by the patient. These analyses are not intended to provide any binding specifications regarding the diagnosis, prevention or treatment of disease and are non-binding suggestions without obligation. They are merely supportive and complementary to the physician's actual decision regarding diagnosis and therapy. Further diagnostic evaluation by the doctor is recommended and other treatment options are not excluded. The freedom of the treating physician in prescribing and selecting therapies shall remain unaffected.

Genetic Investigation Requisition

Patient

Surname	
First name	
Date of Birth	
Gender	☐ female ☐ male
Street	
City, Postcode	
Country	
Telephone	
Email	

Doctor / Clinic

Institution	
Doctor	
Street	
City, Postcode	
Email	

Stamp and signature of responsible medical person:

Samples

☐ FFPE tumor sample	☐ FROZEN tumor sample	☐ EDTA blood (1 x 0,5 ml)
other		
Diagnosis		
Tumor type/ICD-10 code		
Tumor purity (%)		Tumor sample location

I hereby authorise Mr / Mrs

(responsible medical person, in accordance with the German Gene Diagnostics Act)

to conduct a genetic investigation, the sampling required for this and to analyze the transcriptome and exome of the sample for the purpose of a molecular examination of the tumor in relation to normal tissue. I will bear the associated costs as listed in the attached appendix "Service Description", regardless of a possible reimbursement provided by my health insurance.

.....
Date, Signature
(Patient)

Declaration of Consent

With my signature I hereby give my informed consent according to German Gene Diagnosis Act (for content and scope see below) to undergo a genetic investigation and the associated required sampling. I have had the opportunity to ask questions and I am aware of my right to revoke. I agree that samples used for the genetic investigation are given to Alacris Theranostics GmbH and in this respect the responsible medical person is released of his or her duty of confidentiality.

I have been informed about the nature, significance, implications and implementation of genetic testing, possible result-based treatment options and the risks involved, according to §9 of the German Gene Diagnosis Act, on the basis of the attached information sheet, which was given to me and explained.

With regards to the storage and use of study samples and analysis results for scientific purposes, additional tests or quality control, in anonymized form for max. 10 years I :

☐ agree ☐ do not agree

With regards to the storage and use of analysis results beyond the time span of 10 years I:

☐ agree ☐ do not agree

In rare cases, medical findings not related to the investigation mandate can become apparent. Regarding these findings I would:

☐ not ☐ in every case like to be informed

I agree that findings from the genetic investigation can be sent to other doctors/persons (as indicated below).

☐ yes ☐ no

1. Name, Address, Email:

2. Name, Address, Email:

3. Name, Address, Email:

I am aware that my participation in this investigation is voluntary and I may withdraw consent at any time. I acknowledge that a sufficient reflection period was given between being provided with relevant information and providing my consent. I also had the opportunity to ask additional questions.

Date, Signature
(Patient)

Date, Signature
(Responsible medical person)

Dear Patient,

There is a possibility you or a family member may have a hereditary disease, or you are considering prenatal diagnostics or cancer has been diagnosed. Many genetic changes can now be detected by examining an appropriate sample. Before any genetic test is performed your doctor must inform you about its nature, meaning and scope (obligation to inform). In the case of investigations involving cancer patients (spontaneous or hereditary cancer) genetic counselling should be offered; however, you may waive this counselling, in writing, if you have previously been informed about the counselling content. The information below therefore provides an explanation of what happens during genetic counselling and what can be achieved through this process. Genetic counselling is offered to help you to answer questions regarding how a possible hereditary condition or genetic disease impacts you (or your family) and correctly assess potential consequences (risks) for your life and family planning. How far it can actually be useful depends on the particular disorder and your personal questions. Both determine what is discussed in the context of genetic counselling and the accuracy of statements about genetic risks the context of genetic counselling and the accuracy of statements about genetic risks.

Genetic counselling routinely includes:

- Clarifying your personal questions and the objectives of the counselling
- Assessment of your personal and family medical history
- Evaluation of available medical findings and diagnostic reports
- As accurate as possible genetic diagnosis
- Detailed information regarding diseases or disabilities in question
- Assessment of specific genetic risks
- Description of the common genetic risks
- A detailed consultation regarding the possible meaning of the information above, where appropriate, for your health and wellbeing, as well as life and family planning.

Genetic counselling sometimes includes:

- A physical examination of you or family members
- Sample taking (e.g., blood, saliva, amniotic fluid, etc.).

If during the course of the consultation it turns out that there are some additional risks so far unknown to you, you would be informed about it. You can, however, determine if and to what extent you want to be informed and whether you want to undertake further medical-genetic examinations. Genetic tests would not be performed without your positive decision to do so.

During the process of genetic testing, abnormalities may be revealed, which according to the current medical knowledge, are not known to cause adverse health effects. In such cases you will only be notified about these abnormalities if deemed a necessary part of the requested investigation.

A comprehensive overview and diagnosis of all aberrations is not possible, as often the cause is not known. Nor is it possible to exclude any risk of illness for yourself or your family, including your children. In some cases, no exact guidance on the probability of occurrence of a particular disease or disability is possible; even if only a low risk of recurrence is given, this means that an occurrence is possible.

The most important points arising from the consultation will be documented and summarized for you in writing. You decide which doctors are informed. Your blood/tissue samples may be analysed in different laboratories. In such cases everyone involved is subject to all legal requirements of doctor-patient confidentiality and, in particular, data protection compliance.

For further information please do not hesitate to contact us.

Sample Requirements & Shipping Guidelines:

blood and tissue samples



Please note: For successful molecular characterisation of tumors, high quality nucleic acids are essential. This can only be ensured by proper handling and shipment of samples. Please read the guidelines carefully and send samples only after consultation with an Alacris representative.

Sample types required:

1. Tumor tissue for exome and transcriptome sequencing (**with at least 30% tumor purity**)
2. Blood (or healthy tissue) sample for healthy (control) exome sequencing.

The best material for molecular characterisation of tumors is fresh tumor tissue, isolated during OP, immediately frozen on dry ice and shipped to Alacris on dry ice. If there is no fresh tissue available, we can also use FFPE tumor material for analysis.

General note for shipping samples:

Before shipping samples, please contact Alacris to clarify date and type of shipment. If necessary, Alacris can also organize delivery of samples via a courier service. For dry ice shipment, inform the courier two days prior to shipment date.

e-mail: cmtanalysis@alacris.de **Telefon:** +49 30 8431 225 - 10 / - 30

Sample transport containers:

Alacris can supply (free-of-charge) the appropriate consumables required for blood collection, sample transport and storage. If self-supplied consumables are used, please ensure that sample tubes can withstand snap-freezing in liquid nitrogen/dry ice.

Important notes about samples

Fresh tumor sample

Equipment: - 80 °C freezer (for sample storage before shipment), liquid nitrogen, Styrofoam box (or equivalent), dry ice

- Optimal tissue size: (2 mm x 2 mm). Minimum tumor purity is 30 %.
- Transfer the tissue to a cryotube and immediately freeze on dry ice; when immediate freezing is not possible, the sample can be kept on ice for a maximum of 5 min.
- Store the tubes at -80 °C prior to shipping.
- For shipment, put tubes in in a sealed plastic bag, and place in a Styrofoam box (or equivalent) with dry ice (at least 10 cm of dry ice is required above and below the sample tubes).

Important: The completed Tumor Sample(s) Description form should be sent along with the sample.

→ Sheet 2

Sample Requirements & Shipping Guidelines:

blood and tissue samples

FFPE tumor sample

- Optimal size of the FFPE-embedded tissue is 2 cm² × 3 mm with a **tumor content of at least 30%**
- 8–10 unstained tissue sections, each 10 µm thick.
- Sections should be provided either in 1.5 mL tubes, with two sections per tube, **or** unstained on glass slides
- Before collection of the sections, the first 2-3 sections should be trimmed and discarded from the block
- Slide with a stained section of the FFPE block corresponding to the tumor or a photo of the stained section.
- Shipment in a padded envelope with protective casing for the glass slides (if applicable)

Important: The completed Tumor Sample(s) Description form should be sent along with the sample.

Blood sample

We need one tube with 0,5 ml blood. Use only EDTA as anticoagulation agent. Samples with heparin cannot be processed.

- Fill EDTA tube with 0,5 ml blood and label the tube with unique patient number.
- Let the tube stand for 30 min at room temperature before proceeding to Blood shipment variant 1 or Blood shipment variant 2.

Blood shipment variant 1: frozen

Equipment: -80 °C freezer (for sample storage before shipment), Styrofoam box, dry ice

- Freeze the EDTA/blood tubes on dry ice.
- Store the tubes at -80°C prior to shipping.
- For shipment, put tubes in a sealed plastic bag and place in a Styrofoam box (or equivalent) with dry ice (at least 10 cm of dry ice is required above and below the sample tubes).

Blood shipment variant 2: refrigerated

Equipment: Styrofoam box (or equivalent), +4°C cool blocks

- Put the EDTA/blood tubes in a sealed plastic bag and place in a Styrofoam foam box with +4°C cool blocks.
- Send immediately, DNA should be isolated within 24 hours.
- Select a “handle with care” shipment option, shaking should be avoided.

Important: Inform Alacris at least two days before you plan to ship the non-frozen blood.

Please note: Alacris does not carry out pathological assessment of the provided samples. Samples have to be selected by the treating physician/pathology lab that provides the sample. Any dissection (including macro dissection) of the tumour material - including the decision to perform it - must be carried out by the providing pathologist or pathology laboratory.

IMPORTANT:

Analyses cannot start without the pathologist's sample description information below which is mandatory.

Tumor sample(s) description

Please fill in and send to Alacris along with the sample(s) for analysis.

Patient

Name

Date of Birth

Gender ☐ female ☐ male

Contact person

Institution

Name

Telephone

Email

Tumor sample(s)

Sample(s) obtained on [DD/MM/YYYY]

Pathology report(s) ☐ enclosed ☐ will be sent by [DD/MM/YYYY]

For FFPE sample(s) only: ☐ slide(s) with stained tissue sections or photo thereof is included.

Nr.	Sample ID (e.g. FFPE block number)	Tissue type	Tumor content * in [%]	Type of sample e.g. FFPE, frozen, etc.de	Comments

* estimation of the tumor content in the sample is essential for the analysis

For multiple tumor samples only:

The samples were obtained:

☐ from different parts of the same tumor region ☐ from different independent tumor regions

Comments:

Important: Alacris does not carry out pathological assessment of the provided samples. Samples have to be selected by the treating physician/pathology lab that provides the sample. Any dissection (including macro dissection) of the tumour material - including the decision to perform it - must be carried out by the providing pathologist or pathology laboratory.

Checklist of the documentation

required by Alacris

Please note: Processing of sample(s) cannot be started until all listed forms and information are supplied to Alacris.

The patient was given the 'Patient Information Sheet: Genetic Counselling before genetic investigations in accordance with the Gene Diagnostics Act - GenDG, and duly informed

☐

Order / Contract

'Genetic Investigation Requisition' and 'Declaration of Consent' to genetic investigation in accordance with the German Gene Diagnostics Act (form attached) signed by

the doctor ☐ and the patient ☐

☐

Patient information

Information on disease and treatment history

☐

Information about previous genetic tests performed on the patient and their results

☐

Sample information

Pathology report(s) on the sample(s)

☐

Tumor Sample Description (form attached)

☐

Documentation is complete

☐